


## Comparing h2bc and 1,1-adequate

NM160003E

Both h2bc and 1,1-adequate are designed to separate  $^2J_{CH}$  and  $^3J_{CH}$  when a combination of cosy and hmbc is unable to identify correlation signals. Features and best applications of these two techniques, which are used for the same objective, are described below.



4-methylumbelliferone


h2bc can only acquire limited correlation signals from 4-methylumbelliferone rich in quaternary carbon shown in the figure. 1,1-adequate can detect more correlation signals although the process takes time. Also, because the fatty series and aromatic ring have different  $^1J_{CC}$  values, setting measuring conditions using either value will change the levels of sensitivity of the correlation signals detected.

### 【h2bc】


- Correlation of  $2J_{CH}$  only observed
- Higher sensitivity than 1,1-adequate
- No correlation with quaternary carbon observed
- Good result acquired under default conditions

### 【1,1-adequate】

- Correlation of both  $1J_{CH}$  and  $2J_{CH}$  observed
- Lower sensitivity than h2bc
- Correlation with quaternary carbon observed as well
- Correlation signal may not be visible under certain conditions ( $1J_{CC}$  value)



h2bc, 4 scans (25 min)



1,1-adequate, 48 scans (7h)

Instrument: JNM-ECZ400S+ROYAL probe; sample: 100 mg 4-methylumbelliferone