

## Study on phase transition of crystalline polymers by using new extended VT HXMAS probe

NM160009E


Almost all crystalline polymers are mixture of crystalline and amorphous phases. Molecular mobility of these phases which greatly influence the mechanical property of materials changes at a glass transition temperature ( $T_g$ ) and melting temperature ( $T_m$ ). The new extended VT HXMAS probe is broadly applicable to study on thermal properties of crystalline polymers.

The basic methods of solid state NMR, CPMAS and DDMAS, are useful for identification of crystalline and amorphous peaks and analysis of their molecular dynamics.

Here we present practical temperature dependent analyses for crystalline polymers.


### Ex.1) poly(L-lactic acid)

JNM-ECZ600R  
3.2mmHXMAS/VT  
MAS : 5kHz  
Temp. correlation:  $Pb(NO_3)_2$



- \* An additional option may be need to achieve 200°C sample temperature.
- \* A heat-resistance zirconia sample tube is required.

## Ex. 2) poly( $\epsilon$ -caprolactone)



\* A liquid nitrogen dewar vessel is required for the measurements below room temperature.

For the DDMAS measurements, signals from mobile regions are selectively enhanced by setting relatively short relaxation delay, 2s.

Below  $T_g$ , signals are observed only in the CPMAS spectrum because both crystalline and amorphous regions have low molecular mobility. Above  $T_g$ , rubber state amorphous signals appeared in the DDMAS spectra. Above  $T_m$ , very narrow molten state signals are observed in the DDMAS spectrum while any signals are absent in the CPMAS spectrum.