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
➤ 打开Spin

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The screenshot displays the Bruker console software interface. At the top, there are tabs for 'Samples / Jobs / Queue', 'Monitor', and 'Status'. Below this, a status bar shows 'Lock & Shim', 'Sawtooth', and 'Probe' options. The main interface is divided into several panels:

- Sample State:** Includes a solvent input field set to 'NONE' and a lock meter showing a value of 864.
- Spinner:** Features a 'Spinner' icon, a 'Current' field set to 16[Hz], and a 'Target' field set to 15[Hz]. A 'SPIN ON 15[Hz]' indicator is visible.
- Temperature:** Shows a 'Current' field set to 24.9[dC] and a 'Target' field set to 25.0[dC]. A 'TEMP OFF 25[dC]' indicator is present.
- Lock Control:** Contains a 'Status' field set to 'LOCK OFF' and a 'LOCK ON' button. It also includes a checkbox for 'Automatic Gradient Shimming' and a button labeled 'Open the Gradient Shim tool' with a shim icon circled in red.
- Shim Groups:** A section for adjusting shim groups Z1, Z2, Z3, and Z4, each with a '0[Hz]' target and adjustment buttons.

At the bottom of the interface, a status bar displays various parameters: Receiver Gain: 50, Spin: 15[Hz], Lock: 903, Temp: 25[dC], Helium: 50[%], Nitrogen: 75[%], and No Jobs.

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


The screenshot displays the 'Spectrometer Control - Smart Mode' interface. The main window shows a table with columns for 'No.', 'Sample Name', 'Solvent', and 'Comment'. A single row is visible with '1', 'test', and 'NONE'. The right-hand panel, titled 'Job', contains 'Open Jobs' and 'Method Parameters' sections. The 'Method Parameters' section is expanded to show the following settings:

Parameter	Value
wet_solvent	Chloroform
peaks	1
force_tune	<input type="checkbox"/>
autogain	<input checked="" type="checkbox"/>
rgain	<input type="checkbox"/>
scans	8

Below the parameters, there are fields for 'Scans' (set to 8), 'Policy' (Choose a scheduling policy), and 'Start' ([dd-mmm-yyyy] hh:mm:ss). A 'Submit' button is located at the bottom right of the job configuration panel. The bottom status bar shows system information: Receiver Gain: 50, Spin: 0[Hz], Lock: 431, Temp: 25[dC], Helium: 50[%], Nitrogen: 75[%], and Queue Length: 0.

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