

Protocol: Observation of Bacteria Flagellum (Gold Labeling)

Required reagents

- 1) PBS (phosphate buffered saline)
- 2) 4% paraformaldehyde/1% glutaraldehyde in PBS
- 3) 0.01% Poly-L-lysine solution
- 4) DDW (double distilled water)
- 5) Positively charged nanogold (Nanoprobes, #2022)
- 6) Gold Enhance EM Formulation (Nanoprobes, #2113)
- 7) Dextrose in Water (10 mg/ml)


ASEM dish film coating

1. Place 0.01% Poly-L-lysine solution in the ASEM dish. Leave for 30 min at room temperature.
2. Wash three times with DDW.

Sample preparation

1. Centrifuge 1000 μ l of the bacteria suspension at 3000 rpm for 5 minutes. Remove the supernatant and add 1000 μ l PBS. (Repeat three times.)
2. Remove the supernatant and fix with 200 μ l 4% paraformaldehyde/1% glutaraldehyde for 10 minutes.
3. Centrifuge at 3000 rpm for 5 minutes. Remove the supernatant and add 1000 μ l DDW. (Repeat three times.)
4. Remove the supernatant. Suspend bacterial cells in 100 μ l DDW.
5. Apply 20 μ l bacteria suspension to the ASEM dish membrane and incubate for 30 min at room temperature.
6. Wash three times with DDW.
7. Apply 3 μ M positively charged Nanogold in DDW and incubate 20 min.
8. DDW wash; 5 min x 3.
9. Gold enhancement: Nanoprobes, #2113 protocol.* Development time: 10 minutes.
10. Exchange the supernatant for dextrose (10 mg/ml) prior to ASEM observation.
11. Recommended ClairScope conditions: spot size = 20 - 35, acceleration voltage = 30 kV (Fig. 1).

*Nanoprobes #2113 protocol: steps #3 & 4 (PBS wash with 50 mM glycine and PBS-BSA wash) were not used.



x20,000 scale bar = 1 μ m

x20,000 scale bar = 1 μ m

Fig. 1. ASEM images of *Salmonella enterica* NBRC 13245

The bacteria specimen was kindly supplied by Dr. Seiichi Haga, Meijyo University, Japan